



Leak testing services

Individual services in leak testing technology for your applications

Leak testing services

Our services at a glance

Feasibility / cycle-time studies

In the course of a feasibility study we gain an understanding of the technical requirements of your product. This enables us to develop and apply a leak test procedure in industrial tightness control. This applies for new products prior to the market introduction or existing products with the need for cycle-time optimization. We work out an individual optimum solution for you to meet tightness requirements and respective regulations in detail.

Feasibility study for Container Closure Integrity Testing (CCIT)

Quality and effectiveness of drugs depends significantly on their proper packaging to maintain sterility throughout their lifetime. With a feasibility study we gain an understanding of your packaging and the CCIT capabilities. With our three CCIT-technologies we will confirm achievable detection limit and cycle time on your packaging.

Contract leak testing

When investment in leak testing instruments is not an economically meaningful expense, we are your choice by offering contract leak testing as a service.

Residual gas analysis

Residual gas analysis in vacuum systems and tightness control with a wide range of tracer gases like air components or electrolytes. The results are transferred to you in a detailed report.



We provide you the following solutions

Pfeiffer Vacuum is providing complete vacuum solutions. We offer products and technical assistance over the entire process. Further, we offer test procedure development and test implementation. In close cooperation with our customers, we create customized solutions for your product. If you have components or final products, which should be leak tested, we can provide you feasibility studies or cycle-time studies as an all-in-one service in our application laboratory.

- Assistance in determining your leak rate requirements and conversion into a test recipe
- Technical support in selecting the optimal and safest leak detection solution for your application
- Wide range of solutions: Leak detection using tracer gas or air with special solutions for hermetically sealed products
- Integrated all-in-one solutions, including consulting and a network of partners for automated systems



Automotive



Pharmaceutical packaging



IP protection class



Refrigeration and air conditioning

We cover **15 different test methods** according to DIN EN 1779 / DIN EN ISO 20485 and we meet the requirements inside the USP 1207 with our 4 technologies.

Air leak testing

Fast and reliable leak testing with air for packaging and electronics as well as for industrial and medical applications

Tracer gas leak testing

Highest sensitivity and fast leak testing for high end applications such as automotive, medical & semiconductor industry

Pfeiffer Vacuum products



E-PDQ



ASM 340

Test methods per technology

D1: Pressure decay test, manually via gauge

D2: Pressure rise test, manually via gauge

D3: Pressure change test, bell pressure change, with overpressure, manually via gauge

D3: Pressure change test, bell pressure change, with vacuum, manually via gauge

D4: Flow measurement, **outside/in;** from higher pressure to vacuum; Mass Extraction / **inside/out;** from higher pressure to atmospheric; Micro-Flow

A1–A3: Helium leak testing, **A1** total vacuum technique, **A2** partial vacuum technique, **A3** local vacuum technique,

B2: Vacuum box

B3: Pressure technique by accumulation

B4: Sniffing test

B5: Pressurization – evacuation test, “bombing test”

B6: Sealed objects by external vacuum technique, integral vacuum test

B7: Carrier gas technique

Hint: **A1, B2, C3, D4...** = Abbreviation for test methods according to DIN EN 1779 / DIN EN ISO 20485

**Multi gas analysis/
Optical emission
spectroscopy**

Integrity testing for highly demanding pharmaceutical packages and advanced sealed devices



AMI 1000

**Quadrupole
mass spectrometry**

Residual gas analysis in vacuum systems and tightness control with a wide range of tracer gases like air components or electrolytes



PrismaPro

B6: Sealed objects by external vacuum technique, integral vacuum test

Detection of e. g. Air, Nitrogen, Argon, Helium, and many other gases

B6: Sealed objects by external vacuum technique, integral vacuum test

Detection of any gases (1–200 μ)

Examples of test methods



Integral test of parts under vacuum in helium



Integral test of helium filled parts in vacuum



Sniffing test



Micro-Flow (over pressure test)

VACUUM SOLUTIONS

FROM A SINGLE SOURCE

Pfeiffer Vacuum stands for innovative and custom vacuum solutions worldwide, technological perfection, competent advice and reliable service.

COMPLETE RANGE

OF PRODUCTS

From a single component to complex systems: We are the only supplier of vacuum technology that provides a complete product portfolio.

COMPETENCE IN

THEORY AND PRACTICE

Benefit from our know-how and our portfolio of training opportunities!

We support you with your plant layout and provide first-class on-site service worldwide.

Are you looking for a perfect vacuum solution? Please contact us:

Pfeiffer Vacuum
T +49 6441 802-0
leak-testing-services@
pfeiffer-vacuum.de

www.pfeiffer-vacuum.com



(PI 0507 PEN | August 2020/PoD)